

### MARC-145 Cell Culture

### Lin Yuan

China Animal Disease Control Center (CADC)
OIE Reference Laboratory for Porcine
Reproductive and Respiratory Syndrome

2015.4.13

### Cell Culture Basics



- \*Porcine alveolar macrophages (PAM) and MARC-145 cells are known to be permissive to PRRSV infection and replication.
- \*The MARC-145 is routinely used for in vitro propagation of PRRSV and for large scale production of PRRSV vaccine strains.

### Cell Culture Basics



#### MARC 145 Cell line


Cell Lines	<b>MARC 145</b>
Organism	Cercopithecus aethiops
	African green monkey
Tissue	kidney
Cell Type	epithelial
<b>Culture Properties</b>	adherent
Virus Susceptibility	PRRSV

### Cell Culture Basics



#### MARC 145 Cell Culture Environment

Cell Lines	<b>MARC 145</b>
Basal Media	DMEM (Dulbecco's Modified Eagle Medium)
Serum	5-10%FBS (Fetal bovine serum)
PH	7.4
$CO_2$	5%
Temperature	37°C
Subculture Ratio	1:3-1:5
<b>Medium Renewal</b>	Every 2 to 3 days.



### Materials Needed

- Tissue-culture treated flasks, plates or dishes
- Complete growth medium
- Water bath 37°C
- Disposable, sterile pipettes
- Incubator with humid atmosphere of 37 °C 5% CO<sub>2</sub>
- Inverted microscope
- 0.25% Trypsin-EDTA
- Automated cell counter

- Put on gown and gloves, turn on cabinet and swab all internal surfaces of the cabinet with 70% alcohol.
- Prepare Complete Growth Medium for the cell line. Medium should be at 37°C.
- Observe cell monolayer.



- **Remove** and discard the cell culture media.
- Wash the cell layer with 1ml trypsin-EDTA solution
- Add 1ml trypsin-EDTA solution to the flask
- **Incubate** the culture vessel at 37°C for 5-15 minutes




- **Observe** the cells by inverted microscope for detachment.
- **Inactivate** trypsin with 4ml complete growth medium





- Count the cells in suspension
- Inoculate cells to culture vessals



## Cell Growth and Propagations Company C

- **Incubate** the cells at 37°C 5% CO<sub>2</sub>
- Every 2 to 3 days, subculture ratio is 1:3 to 1:5.



# Cell Growth and Propagations Company of the Company

### Precautions

- Cells should be passaged when they are in the log phase
- Prewash cells is necessary.
- Prewarming of the enzyme solution will decrease the exposure period
- Trypsin need to inactivate with serum.
- To be sure that most (>95%) are single cells.
- Always use proper aseptic technique and work in cell culture hood.
- Always wear personal protective equipment

