

Antimicrobial susceptibility testing (Broth microdilution method)



Rakuno Gakuen University
School of Veterinary Medicine



WOAH Collaborating Centre
for Food Safety




World Organisation
for Animal Health
Founded as OIE




2025. Jul.

**Rakuno Gakuen University
Veterinary AMR Center, NVAL, MAFF**

Antimicrobial susceptibility test method




Reagents and equipment prepared by the broth microdilution method



- ① Microtiter plate
- ② Lid
- ③ Reading mirror
- ④ Tip
- ⑤ 12 or 8 channel
Pipette
- ⑥ Micropipette
- ⑦ Tested strains
- ⑧ Reservoir
- ⑨ Cation adjusted
Muller-Hinton Broth
- ⑩ 3 mL sterile saline
- ⑪ Turbidity standard
solution

Adjustment turbidity of bacterial suspension




- Pick a test organism on an overnight agar culture
- With 3 mL of sterile saline suspension (Adjust the turbidity so that it becomes turbid between Mcfarland standard turbidity solution 0.5 and 1.)

Adjustment of inoculating bacterial suspension



**Add 0.025 mL (25 μ L) of the previously prepared bacterial suspension in addition to 12mL Muller-Hinton broth
Mix evenly and use this as the bacterial suspension for inoculation.**

Inoculation and culture of the bacterial suspension



- Transfer inoculum into reservoir
- Inoculate the wells of the microtiter plate with 100 μ L of the bacterial suspension using an 8- or 12-channel pipette fitted with tips, starting with the wells with the lowest drug concentration.

✖ Remove the microtiter plate from the refrigerator beforehand and bring it back to room temperature.

Since the antimicrobials is dry and stuck to the bottom of the wells of the dry plate, be careful to avoid touch with the tips during inoculation.

- After inoculation, cover the plates and incubate under the culture conditions specified for each bacterial species. (When culturing, the number of stacked plates should be limited to a maximum of three.)



Evaluation


[Positive Growth Criteria]

- Turbidity or precipitation of 1 mm or more in diameter is observed
- Two or more precipitates are found even if the diameter of the precipitates is less than 1mm


[Negative Growth Criteria]

- No visible turbidity or precipitation
- 1 precipitate less than 1 mm in diameter, even if present


Positive



Positive




Positive




Turbidity

Negative



Judgment example



Minimum inhibitory concentration (MIC)

	Drug dilution concentration ($\mu\text{g/mL}$)								
	0.5	1	2	4	8	16	32	64	Control
Growth of bacteria	+	+	+	+	+	-	-	-	+

MIC 16 $\mu\text{g/mL}$

	Drug dilution concentration ($\mu\text{g/mL}$)								
	0.5	1	2	4	8	16	32	64	Control
Growth of bacteria	+	+	+	+	+	+	+	+	+

MIC > 64 $\mu\text{g/mL}$

	Drug dilution concentration ($\mu\text{g/mL}$)								
	0.5	1	2	4	8	16	32	64	Control
Growth of bacteria	-	-	-	-	-	-	-	-	-

MIC $\leq 0.5 \mu\text{g/mL}$

Re-test required

Drug dilution concentration	Drug dilution concentration ($\mu\text{g/mL}$)								
	0.5	1	2	4	8	16	32	64	Control
Growth of bacteria	+	+	+	+	+	-	-	-	-

In the absence of control growth

Drug dilution concentration	Drug dilution concentration ($\mu\text{g/mL}$)								
	0.5	1	2	4	8	16	32	64	Control
Growth of bacteria	+	-	+	-	+	-	-	-	+

Growing (skipping) in excess of concentration

For QC strains

- When the micro dilution method is performed, the MICs of the **QC strain (according to CLSI)** are also measured.
- The MIC of the QC strain must be within the range of the QC reference value.